

\\USER\MPI\praveen\20230119_TRTE\pivGRE_B0map_R2_1mm

TA: 0:41 PM: FIX Voxel size: 1.1×1.1×1.0 mmPAT: 2 Rel. SNR: 1.00 : 47791c8

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R9.4 P1.6 H20.5 mm
Orientation	T > C-41.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	20.0 %
Slices per slab	30
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	10.0 ms
TE 1	2.27 ms
TE 2	4.54 ms
TE 3	7.11 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	R01-31

Contrast - Common

TR	10.0 ms
TE 1	2.27 ms
TE 2	4.54 ms
TE 3	7.11 ms
MTC	Off
Magn. preparation	None
Flip angle	8 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution - Common

FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
Base resolution	200

Resolution - Common

Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Accel. factor 3D	1
Reference scan mode	Integrated

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R9.4 P1.6 H20.5 mm
Orientation	T > C-41.4
Phase enc. dir.	A >> P
Slice oversampling	20.0 %
Slices per slab	30
FoV read	210 mm
FoV phase	100.0 %
Slice thickness	1.00 mm
TR	10.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	R9.4 P1.6 H20.5 mm
Orientation	T > C-41.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R9.4 P1.6 H20.5
R	9.4 mm
P	1.6 mm
H	20.5 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-41.4
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None

Geometry - Saturation

Special sat.	None
--------------	------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

Geometry - Tim CT

Tim CT mode	Off
Slabs	1
Slices per slab	30
Slice thickness	1.00 mm
Dist. factor	20 %
FoV read	210 mm
FoV phase	100.0 %
Segments	1

System - Miscellaneous

Positioning mode	FIX
Table position	F
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R6.5 P19.5 H6.5 mm
! Orientation	T > C-40.9
! Rotation	0.00 deg
! A >> P	112 mm
! R >> L	200 mm
! F >> H	23 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	399.717497 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
-----------------	------

Physio - Signal1

TR	10.0 ms
Concatenations	1
Segments	1

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	210 mm
FoV phase	100.0 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	1

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	1
Contrasts	3
TR	10.0 ms
TE 1	2.27 ms
TE 2	4.54 ms
TE 3	7.11 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	Off
Phase stabilisation	Off
Asymmetric echo	Off
Contrasts	3
Flow comp. 1	No
Readout mode	Monopolar

Sequence - Part 1

Multi-slice mode	Interleaved
Bandwidth 1	890 Hz/Px
Bandwidth 2	890 Hz/Px
Bandwidth 3	890 Hz/Px

Sequence - Part 2

Segments	1
Acoustic noise reduction	None
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On

Sequence - Special

RF Pulse Duration	1000 us
RF Pulse TBW	9.0
Invert Slice Sel. Grad.	On

Sequence - Assistant

Mode	Off
------	-----